

# ERGO-Lux MP mobil

The new star in the sky of magnifiers



**Our most versatile  
series of illuminated  
hand magnifiers!**

# Innovations made by SCHWEIZER

SCHWEIZER, innovation leader for illuminated magnifiers, has provided solutions for people with impaired vision for decades.

The low vision specialist was the first company worldwide to develop a modular magnifier system. Additional world premieres included the first magnifiers with LED illumination and, above all, the first comprehensive ergonomic magnifier series. All the new products represented milestones in low vision care and have meanwhile become proven aids for improving visual acuity.

## ■ MODULAR

First modular magnifier system



2000

## ■ ÖKOLUX mobil

First LED illuminated hand magnifier



2004

1990

2000

## ■ ÖKOLUX

First LED illuminated stand magnifier



## ■ MODULAR AS

Illuminated magnifier with automatic sensor



# Our mission is to aid people who are visually impaired!

As a result of the changing age structure of our population and hence gradual loss of sight, we need to continuously adjust our range of magnifiers and other visual aids. Older people have increasing difficulty with contrast sensitivity and movement coordination. In the end, they require magnifying visual aids that help them master everyday tasks.

Since the causes and effects of visual impairment vary widely, the solutions offered must be tailored to the individual needs of those with low vision: a wide range of power variations and different fields of view to meet individual requirements should be available.

The size of the visual aid is also key to success when it comes to different applications. Depending on the customer, the design can play a major role in the decision to purchase...



Alfred Schüttinger, Managing Partner of A. SCHWEIZER GmbH

## ■ ERGO-Line

First ergonomic range of magnifiers



2009

## ■ ERGO-Lux MP

First illuminated magnifier with boost switch



2015

2007

2013

## ■ ÖKOLUX plus

First series of magnifiers in 3 light temperatures



## ■ ERGO-Lux MP mobil

The new star in the sky of magnifiers





# New: ERGO-Lux MP mobil – small yet the perfect complement!

The new ergonomic illuminated hand magnifier combines many special features which make it an indispensable tool for visually impaired people. A true bestseller thanks to its innovative design.

## ■ Versatile solutions

- 14 magnifiers
- 11 power variations
- 3 light temperatures (warm white, cool white, neutral white)
- 2 brightness levels

With this series of portable, ergonomic magnifiers, eye care professionals can serve a wide range of customers. In addition to the low vision market and people wearing glasses, they can also reach those who want to benefit from a magnifier despite their usually good vision.



# ERGO-Lux MP mobil – handy & extra bright


## ■ Maximum comfort thanks to design and material

Featuring a handle that fits comfortably in the hand since its shape promotes a natural wrist posture and its prominent handle end supports the ball of the wrist while ensuring reliable scrolling on the object. Thumb and index finger coordinate the fine movements and automatically find the thumb and index finger rests. Extra grip strength is added by an anti-slip material consisting of 2 different textures — a pleasant soft-touch zone.

## ■ Carefully designed control elements

- Battery-operated LEDs provide superior, shadow-free lighting
- **Multi-Power** light option ensures adjustment to individual lighting requirements
- Slightly raised slide-switch for user-friendly, haptic guidance
- Switch operation provides audible feedback
- Battery box closes with distinct click for quick and easy battery replacement

- In boost mode: Brightest, effective lighting
- Unique, long-lasting consistent light intensity – no loss of brightness (see chart)



## ■ Switch position 1: regular illumination

- Energy-saving
- Homogeneous, even spread of light
- For normal ambient lighting

**Our brightest magnifier!**

## ■ Switch position 2: boost mode

- 50% greater brightness compared to standard lighting
- Excellent contrast for difficult viewing tasks and low-contrast objects
- Particularly suitable in the event of insufficient ambient lighting

**ERGO-Lux MP mobil – the stylish, small illuminated hand magnifier with optimised light adjustment for brightness beyond regular lighting...**





# ERGO-Line – The complete range at a glance

## ■ ERGO-Line – the lifestyle magnifier range with the extra edge

Application and aesthetics as an experience: ergonomic design, clear appearance and pleasant feel.

- Form and function perfectly attuned
- Thumb and index finger automatically control the fine movements
- Ergonomic shape in timeless white/charcoal grey
- Suitable for left and right-handed use
- Aspheric lightweight lenses

## ■ ERGO-Lux MP mobil

Portable lifestyle magnifier with **Multi-Power** light option: one of the brightest magnifiers on the market offering 2 brightness levels and 3 light temperatures.


- Boost switch provides 2 brightness levels for optimised lighting
- One of the brightest magnifiers on the market!
- 3 individually tested light temperatures improve reading performance (warm white 2700K, neutral white 4500K, cool white 6500K)
- Low energy consumption and easy battery replacement thanks to simple closing mechanism in the handle
- 14 different lenses with 11 powers



■ ERGO-Lux MP mobil



■ ERGO-Lux MP with ERGO-Base



## ■ ERGO-Lux MP

Functional lifestyle magnifier with **Multi-Power** light option: available in a choice of 3 light temperatures with 2 brightness levels.

- Boost switch provides 2 brightness levels for optimised lighting
- 3 individually tested light temperatures improve reading performance (warm white 2700K, neutral white 4500K, cool white 6500K)
- Low energy consumption and easy battery replacement thanks to simple closing mechanism in the handle
- 5 different lenses with 5 powers

## ■ ERGO-Base

The ERGO-Base fixed-focus holder converts the ERGO-Lux MP into a perfect illuminated stand magnifier.

- Patented magnetic design for easy connection without wear and tear!
- Easy scrolling on the object
- Comfortable reading position thanks to tilt feature of the holder
- Available for all illuminated magnifiers of the ERGO-Lux MP series

## ■ ERGO-Style

Traditional reading magnifier without illumination as functional lifestyle magnifier.

- Ergonomic design and correct distance between the object and the lens ensure a natural and stress-free position
- 5 different lenses with 5 powers

## ■ ERGO-Pocket

Easy-to-use pocket magnifier in an appealing design.

- Fits nicely in any handbag, shirt or trouser pocket
- Convenient integrated lanyard loop
- A single movement of the hand reveals the lens
- 1 magnifier with 6 different power variations



■ ERGO-Style



■ ERGO-Pocket

# ERGO-Lux MP mobil – at a glance



Portable lifestyle magnifier with **Multi-Power** boost option. Outstanding **Multi-Power** benefits: boost switch for superior illumination, choice of 3 light temperatures (warm white – 2700K, neutral white – 4500K and cool white – 6500K).

- Easy-to-feel boost switch provides audible feedback when operated and controls 2 brightness levels
- High contrast, even spread of light thanks to ideal adjustment of light intensity
- Individually tested light temperatures improve reading performance
- Low energy consumption and easy battery replacement thanks to simple closing mechanism in the handle which locks with a distinct click
- Ergonomic design and soft-touch zone for extra grip strength and pleasant feel in timeless white/charcoal grey
- 14 different lenses with 11 powers, 3 light temperatures and 2 brightness levels
- Aspheric lightweight lenses



Even distribution of light in any situation – choice of 3 light temperatures.



Extra bright lighting – the boost switch controls 2 brightness levels.



Easy battery replacement: open battery box by folding back the bottom shell of the handle.



View finest details in a comfortable reading posture.

15-002/3113758/00