

# ERGO-Lux *i* mobil

Ultimate convenience by using  
wireless power



**World  
premiere!**

Light temperatures: 2700 K 4500 K 6500 K



152562 154562 156562  
56 D / Ø 35 mm (opt. Ø 27.0 mm)

152482 154482 156482  
48 D / Ø 35 mm (opt. Ø 27.5 mm)

152392 154392 156392  
39 D / Ø 35 mm

152282 154282 156282  
28 D / Ø 35 mm



152242 154242 156242  
24 D / Ø 55 mm (opt. Ø 52.5 mm)



152202 154202 156202  
20 D / Ø 55 mm



152172 154172 156172  
16 D / Ø 35 mm



152162 154162 156162  
16 D / Ø 60 mm



152132 154132 156132  
12 D / Ø 35 mm



152122 154122 156122  
12 D / Ø 70 mm

# LED illuminated hand magnifiers

**NEW**

## ■ ERGO-Lux i mobil

ERGO-Lux i mobil – the lifestyle magnifier with inductive charging. In addition to 2 brightness levels and a choice of 3 light temperatures, the key advantage of this new magnifier lies in its wireless charging capability which eliminates the need to replace batteries.

- Power supply by induction – no battery replacement, no cable inserted into the magnifier
- Easy-to-feel boost switch provides audible feedback when operated and controls 2 brightness levels for optimised lighting
- High contrast, even spread of light thanks to ideal adjustment of light intensity
- Individually tested light temperatures improve reading performance
- Low energy consumption
- Ergonomic design and soft-touch zone for extra grip strength and pleasant feel
- 17 different lenses with 11 powers, 3 light temperatures
- Rectangular magnifier heads for left- and right-handed use
- Aspheric lightweight lenses
- Order any of these magnifiers as a complete set including charging pad, USB cable and USB power adapter by adding -1 to the item no. (e.g. 152122-1)





Light temperatures: **2700 K** **4500 K** **6500 K**



Box ERGO-Lux i mobil with 2 pads and USB adapters

**150202** **150402** **150602**

Contents: 16 D Ø 60 mm, 10 D 75x50 mm right-handed, 8 D Ø 85 mm, 2x charging pad with USB cable and USB power adapter



**150442** Universal display with ERGO-Lux i mobil

Black coated metal display with exchangeable A4-sized insert and 3 magnetic holders; approx. size (WxHxD): 30.0x30.0x21.5 cm  
Contents: 16 D Ø 60 mm, 10 D 75x50 mm right-handed, 8 D Ø 85 mm (all in 4500 K), 2x charging pad with USB cable and USB power adapter

**164052** Universal display

Black coated metal display with exchangeable A4-sized insert and 3 magnetic holders; approx. size (WxHxD): 30.0x30.0x21.5 cm



**150002** Charging pad

Approx. 70 x 70 mm with 1.5 m USB cable

**150012** USB power adapter

for ERGO-Lux i mobil charging pad

**180902** Pack of 5 lanyards

**180912** Pack of 5 wrist straps



**152112** **154112** **156112**

10 D / 75x50 mm, straight



**152102** **154102** **156102**

10 D / 75x50 mm, right-handed



**152142** **154142** **156142**

10 D / 75x50 mm, left-handed



**152092** **154092** **156092**

8 D / Ø 85 mm



**152082** **154082** **156082**

8 D / 100x75 mm, right-handed



**152072** **154072** **156072**

8 D / 100x75 mm, left-handed



**152062** **154062** **156062**

6 D / Ø 100 mm

# Technology to enjoy...

Maximum user comfort!

## ■ Setting new standards

Innovations made by SCHWEIZER have always been milestones in low vision care and set the pace in the industry.

## ■ The magic word is “induction”

This technology is becoming increasingly popular in everyday devices like smartphones, tablets or wearables.

## ■ Induction charging - how it works

The device is simply placed on the charging pad and then, as if by magic, the integrated battery is powered up without any wire connection between charging station and device. The energy is transferred by induction.

**SCHWEIZER has selected the Qi standard since it is compatible with many devices (mobile phones and other small appliances) that can use the same charging station.**

## ■ Key benefits of ERGO-Lux i mobil

### ■ No more changing of batteries

Cut the cord and take advantage of this convenient feature

### ■ At home:

**Plug in the charging pad once and for all – that’s it!**

The magnifier has its fixed place and is always ready for use – it couldn’t be easier!

### ■ On the go:

**Charging pad with USB port**

Multiple power sources can be used (car, laptop etc.)

### ■ One for all, all for one...

One charging pad suitable for powering up different devices

**Goodbye battery anxiety. Benefit from this innovative technology that eliminates the hassle of changing batteries and inserting a charging cable into the magnifier.**

