

High addition spectacles

for monocular
or binocular
close distance
magnification



OUT NOW
from ImproVision

ImproVision
LowVision Solutions

Prismatic half-eye spectacles

■ IV PRISMATIC

Prismatic half-eye spectacles with a pair of high addition lenses including convergence supporting prisms for binocular close distance magnification.

- Spherical CR 39 lenses
- Black matt PA polymer frame in size 43-25
- Frame size and power printed in temple
- Supplied with protective case
- Made in the European Union



IV PRISMATIC

IV80200 +3.00 D	IV80203 +6.00 D
IV80201 +4.00 D	IV80204 +7.00 D
IV80202 +5.00 D	IV80205 +8.00 D

IVBOXHA Trial box IV PRISMATIC/IV MICROSCOPIC

Content: 1 each of IV PRISMATIC/IV MICROSCOPIC, total 12 half-eye spectacles

IVBOX12 Trial box IV High Additions

Content: 12 IV PRISMATIC/IV MICROSCOPIC spectacles of your choice



Microscopic half-eye spectacles

■ IV MICROSCOPIC

Microscopic half-eye spectacles with one high addition lenticular lens and one CR 39 plano lens on the other side for monocular close distance magnification.

- Spherical CR 39 lenses
- Black matt PA polymer frame in size 43-25
- Frame size and power printed in temple
- Supplied with protective case
- Made in the European Union



IV MICROSCOPIC

IV80220	+ 12.00 D OD	IV80221	+ 12.00 D OS
IV80222	+ 16.00 D OD	IV80223	+ 16.00 D OS
IV80224	+ 20.00 D OD	IV80225	+ 20.00 D OS

IVBOXHA Trial box IV PRISMATIC/IV MICROSCOPIC

Content: 1 each of IV PRISMATIC/IV MICROSCOPIC, total 12 half-eye spectacles

IVBOX12 Trial box IV High Additions

Content: 12 IV PRISMATIC/IV MICROSCOPIC spectacles of your choice



Do you have questions ...

... about our products or do you need further information?
Please give us a call, we have the right specialists for you!

Phone

+49-9191-72 10 56

Fax

+49-9191-72 10 73

E-mail

info@improvision-lvs.com

Internet

www.improvision-lvs.com

15-009/313798/00