

# Technical data – MULTILIGHT Pro

## MULTILIGHT Pro Table Akku 3in1 Rechargeable battery, 3 light temperatures

941103

Foot dimensions (WxDxH): approx. 13.5 x 17.5 x 12.0 cm  
Arm segment length (hinge-to-hinge): approx. 30 cm  
Weight: approx. 1.62 kg; charging cable 150 cm  
All 3 light temperatures individually selectable  
Lithium-ion battery

## MULTILIGHT Pro Table Akku Rechargeable battery

941123 941143 941163

Foot dimensions (WxDxH): approx. 13.5 x 17.5 x 12.0 cm  
Arm segment length (hinge-to-hinge): approx. 30 cm  
Weight: approx. 1.62 kg; charging cable 150 cm  
Lithium-ion battery

## MULTILIGHT Pro Table Mains operated

941023 941043 941063

Foot dimensions (WxDxH): approx. 13.5 x 17.5 x 12.0 cm  
Arm segment length (hinge-to-hinge): approx. 30 cm  
Weight: approx. 1.52 kg; power cord: 150 cm

## MULTILIGHT Pro Table XL Mains operated, extended arm

942023 942043 942063

Foot dimensions (WxDxH): approx. 13.5 x 17.5 x 12.0 cm  
Arm segment length (hinge-to-hinge): approx. 30 cm  
Weight: approx. 3.65 kg; power cord: 150 cm

## MULTILIGHT Pro Table XL Clamp Mains operated, extended arm, clamp

943023 943043 943063

Arm segment length (hinge-to-hinge): approx. 30 cm  
Maximum clamp opening: 5.3 cm  
Weight: approx. 1.2 kg; power cord: 150 cm

## MULTILIGHT Pro Floor 3in1 Mains operated, 3 light temperatures

944003

Lamp pole height: approx. 96 cm; base Ø: approx. 28 cm  
Arm segment length (hinge-to-hinge): approx. 30 cm  
Weight: approx. 3.8 kg; power cord: 150 cm  
All 3 light temperatures individually selectable

## MULTILIGHT Pro Floor Mains operated

944023 944043 944063

Lamp pole height: approx. 96 cm; base Ø: approx. 28 cm  
Arm segment length (hinge-to-hinge): approx. 30 cm  
Weight: approx. 3.8 kg; power cord: 150 cm

2700 K 4500 K 6500 K Light temperatures



# MULTILIGHT Pro

## Light colour in a striking design



20-005 / 370458-EN / 00



# MULTILIGHT Pro

## Best lighting for best vision

The innovative MULTILIGHT Pro low vision lamp series by SCHWEIZER has been developed with the needs of visually impaired people in mind and features numerous impressive details.

Sophisticated technology and maximum durability thanks to selected materials combined with an elegant design and intelligent functionality:  
**MADE IN GERMANY**



Table lamp with rechargeable battery and 3 light temperatures



Table lamp with rechargeable battery in one individual light temperature of choice



Table lamp mains operated in one individual light temperature of choice

## Functionality and design that delight

Designed for visually impaired people, the MULTILIGHT Pro LED lamp series features numerous innovative ideas.

### Design

- 7 different product versions for multiple applications
- Black matte surface prevents irritating reflections
- Robust plastic housing facilitates cleaning
- Elegant and timeless design

### Functionality

- Large central control panel in contrasting colour for easy one-touch operation
- Audible ON/OFF feedback
- Flexible positioning of lamp ensures optimal lighting while minimising glare
- Maximum user comfort thanks to large, easy-to-operate locking levers
- Available in 3 light temperatures (warm white – 2700 K, neutral white – 4500 K and cool white – 6500 K)
- Lamp head can be swivelled in any direction and even pulled right down to the table top
- Recessed connector socket and round barrel connector facilitate plug-in
- Foldable design for space-saving storage; standard arm versions also ideal for mobile use

### Technology

- Low energy consumption and high performance
- Features the tried and tested LEDs of SCHWEIZER's illuminated magnifiers – no mixing of light temperatures
- 2 brightness levels (100% and 70%) for optimal lighting
- Lithium-ion battery with 4-hour operating time and 3-hour charging time
- Easy-to-recognise battery status

