

Prismatic half-eye spectacles

■ BINO med S2

Prismatic half-eye spectacles with a pair of high addition lenses including prisms to support convergence for binocular close distance magnification.

- Spherical CR 39 lenses
- Metal frame in 2 colours in size 40-25, temple length 142 mm
- Supplied with protective case
- Made in Germany



BINO med S2 – silver metal frame

IMP430160	+3.00 D	IMP460160	+6.00 D
IMP440160	+4.00 D	IMP470160	+7.00 D
IMP450160	+5.00 D	IMP480160	+8.00 D



BINO med S2 – gold metal frame

IMP430170	+3.00 D	IMP460170	+6.00 D
IMP440170	+4.00 D	IMP470170	+7.00 D
IMP450170	+5.00 D	IMP480170	+8.00 D