

eMag 70 HD

Comfort in every situation



Technical data

eMag 70 HD

Portable large-screen video magnifier with unique pivotable camera for multiple applications and with a relaxed reading function. eMag 70 HD for extra comfort and flexibility.

- Item No. 571002
- Magnification 2x to 24x, continuous zoom
- Viewing options True full colour, black on white, white on black, yellow on blue, yellow on black + 15 additional contrast modes
- Weight 510 g
- Size Size (WxDxH)
approx. 211 x 147 x 25 mm (folded)
- Operating time Approx. 4 hours
- Battery Built-in rechargeable lithium-ion battery (approx. 3 hours charge time)
- Camera 5 megapixel HD autofocus camera
- Comfort Camera Continuously rotatable camera for near, intermediate and distance vision as well as self viewing
- Display 18 cm (7") LCD display 16:9
- Features
 - Snapshot function with zoom and scroll function, 105 images can be saved and exported via USB
 - Reading line and text masks
 - LED battery indicator
 - LED turn-off option
 - HDMI connection to TV
 - Built-in reading stand (40° tilt)
 - Energy saving mode after 30 minutes
- Accessories
 - Carrying case (WxDxH: 225 x 155 x 45 mm)
 - AC power adapter
 - HDMI cable (2.0 m)
 - USB cable
 - Manual



A. SCHWEIZER GmbH
Hans-Böckler-Str. 7, 91301 Forchheim, Germany
Phone: +49-9191-721056, Fax: +49-9191-721073
info@schweizer-optik.de, www.schweizer-optik.com