

# eMag 50 HD

The HD all-round talent amongst portable video magnifiers



# Technical data

## eMag 50 HD

The portable video magnifier in HD quality has a two-position flip-out handle and an integrated reading stand holding the screen at a natural angle.

- Item No. 550002
- Magnification 2x to 14x continuous zoom
- Viewing options Full colour (with quick change feature), black on white, white on black, yellow on blue, yellow on black  
+ 15 additional contrast modes
- Weight 298 g
- Size Size (WxDxH)  
approx. 139x89x20 mm (closed)
- Operating time Approx. 3 hours
- Battery Built-in rechargeable lithium-ion battery  
(approx. 3.5 hours charge time)
- Camera 5 megapixel HD autofocus camera
- Display 12.7 cm (5") LCD display 16:9
- Features
  - Snapshot function, 80 images can be saved and exported via USB
  - Reading line and text masks
  - LED battery indicator
  - LED illumination can be deactivated
  - Energy saving mode after 30 minutes
- Accessories
  - Case with belt loop and carrying strap
  - Wrist strap
  - AC power adapter
  - USB cable
  - Manual
- Extras Item 544042: Writing and reading stand



A. SCHWEIZER GmbH  
Hans-Böckler-Str. 7, 91301 Forchheim, Germany  
Phone: +49-9191-721056, Fax: +49-9191-721073  
info@schweizer-optik.de, www.schweizer-optik.com