

# eMag 240 HD/200

The classic desktop video magnifiers



14-025 / 313308 / 00

**Comfort  
&  
Clarity**

 **SCHWEIZER**  
Der LowVision-Spezialist

# eMag 240 HD/200

## Technical Specifications

### ■ eMag 200

■ Item No.	520002
■ Magnification	2,3x to 60x
■ Viewing options	Full colour Black on white White on black Yellow on blue Yellow on black + 28 additional contrast modes
■ Weight	18,6 kg
■ Total size	565 x 439 x 480-630 mm
■ XY table size	437 x 406 mm
■ Monitor	20 inch display 16:9 Raises, lowers, tilts and swivels 180 degrees side to side
■ Features	Focus lock Find feature Freeze frame Reading line Illumination can be deactivated Holder for round objects 21 cm working space
■ In the box	eMag 200 AC power supply cable Documentation
■ Accessories	Item 510002: eMag Table high 800 x 600 x 780-1000 mm Item 510012: eMag Table low 800 x 600 x 650-910 mm

### ■ eMag 240 HD

■ Item No.	524002
■ Magnification	1,9x to 71x
■ Viewing options	Full colour Black on white White on black Yellow on blue Yellow on black + 28 additional contrast modes
■ Weight	19,8 kg
■ Total size	565 x 439 x 505-675 mm
■ XY table size	437 x 406 mm
■ Monitor	24 inch display 16:9 Raises, lowers, tilts and swivels 180 degrees side to side
■ Features	Optimized HD quality Focus lock Find feature Position locator Freeze frame Reading line Adjustable text lines & masks Illumination can be deactivated Holder for round objects 21 cm working space
■ In the box	eMag 240 HD AC power supply cable Documentation
■ Accessories	Item 510002: eMag Table high 800 x 600 x 780-1000 mm Item 510012: eMag Table low 800 x 600 x 650-910 mm

