


eMag HD flex

The versatile desktop video magnifier
in high-definition quality!



High contrast full colour

True full colour

Full colour
viewing option with
high contrast

eMag HD flex

Technical Specifications

The eMag HD flex is a self-contained and easy to carry portable desktop video magnifier in a sleek design with an integrated handle. The camera can rotate on all 3 axes to focus on anything around the user.

- Viewing options
 - True full colour
 - High contrast full colour
 - Black on white
 - White on black
 - Yellow on blue
 - Yellow on black
 - + 27 additional contrast modes

- Monitor
 - Adjustable tilt for a comfortable posture
- Power supply
 - 100 – 240 V AC
- Features
 - Autofocus HD camera
 - Focus lock
 - 330 degree camera rotation
 - Distance, document and self view mode
 - Memory of last settings for each mode
 - Find feature
 - Freeze frame
 - Adjustable reading line, mask and shade
 - 41 cm working space

- In the box
 - eMag HD flex
 - AC power supply cable 2.5 m
 - Pocket sized RF remote control
 - Rolling case (not for 24 inch model)
 - Documentation



520022 eMag 200 HD flex

Magnification: 1,4x to 106x
 Weight: 6,2 kg
 Total size: Approx. 466 x 190 x 530 mm
 Monitor: 19,5 inch
 Inclusive rolling case 510 x 280 x 560 mm

522022 eMag 220 HD flex

Magnification: 1,5x to 118x
 Weight: 6,7 kg
 Total size: Approx. 508 x 190 x 530 mm
 Monitor: 21,5 inch
 Inclusive rolling case 510 x 280 x 560 mm

524022 eMag 240 HD flex

Magnification: 1,7x to 131x
 Weight: 7,3 kg
 Total size: Approx. 563 x 190 x 530 mm
 Monitor: 24 inch
 Without rolling case

519022 XY table

Size: 440 x 380 mm

510002 eMag table high

Trapezium shaped with lockable castors, height adjustable; Size: 800 x 600 x 780–1000 mm

510012 eMag table low

Trapezium shaped with lockable castors, height adjustable; Size: 800 x 600 x 650–910 mm