

## ERGO-Line – The complete range


ERGO-Line is an innovative product range that sets standards in terms of form, function and quality. The philosophy is to produce visual aids that focus on the intended use and meet the requirements of visually impaired people while featuring a stylish design.

Enlarging small print texts for better reading, exploring the details of images and identifying structures is what many people want. The ERGO-Line magnifier range ensures all these needs are met while providing a relaxed experience.

The handles of the ergonomic magnifiers combine 4 benefits:

- custom fit in the hand,
- comfortable angle between hand and magnifier,
- additional stabilizing support and
- balanced weight distribution.

All LED illuminated magnifiers of the ERGO-Line range are available in 3 light temperatures:



## Good ideas for improved vision

At home or on the go, for hobbies, leisure and at work: SCHWEIZER magnification and visual aids for multiple applications improve the quality of everyday life.

The range of products includes:

- Magnifiers ■ Video magnifiers ■ High Additions
- Lamps ■ Filters ■ Telescopes

As a leading brand, SCHWEIZER combines experience, expertise as well as innovation. Producing exclusively in Germany, the specialist for magnifiers provides custom solutions for people with impaired vision.

Your partner for improved vision:

## ERGO-Line The ergonomic range of magnifiers ...



... ensures relaxed reading comfort:  
stylish design, excellent magnification  
and superior lighting!

15-001 / 313748/00

## ■ ERGO-Lux MP mobil

The ergonomically designed magnifier handle fits comfortably in the hand. A combination of different materials adds extra grip strength and ensures reliable tracking across the text while providing a pleasant feel. An easy-to-feel, differently coloured boost switch controls the **Multi-Power** light option with 2 brightness levels.

The adjustable LED illumination enables comfortable and relaxed vision even in poor lighting conditions. Lightweight, handy and portable for use in many situations. A reliable aid for hobby and leisure activities as well as in everyday life.

- Boost switch for optimised lighting, with audible feedback
- Choice of 3 light temperatures
- High contrast, even spread of light
- Magnification from 1.5X – 14X (6 D – 56 D)
- 14 different lenses with 11 powers
- Convenient integrated lanyard loop



**Our most versatile magnifier!**



## ■ ERGO-Lux MP

The battery-operated ERGO-Lux MP LED illuminated magnifier with **Multi-Power** light option can always be adjusted to the individual lighting requirements: 50% higher light intensity thanks to the boost switch with 2 brightness levels.

- Boost switch for optimised lighting, with audible feedback
- Choice of 3 light temperatures
- High contrast, even spread of light
- Magnification from 2X – 6X (8 D – 24 D)
- 5 different lenses with 5 powers
- Optional ERGO-Base fixed focus holder

## ■ ERGO-Style

The traditional reading magnifier of our ERGO-Style series comes in a stylish design. Thanks to the ergonomic shape, clear appearance and pleasant feel, its application and aesthetics make a real difference to traditional magnifiers.

- Magnification from 2X – 6X (8 D – 24 D)
- 5 different lenses with 5 powers



## ■ ERGO-Pocket

An attractive accessory for everyone. Highly versatile, it combines practical use with an appealing design. The small yet effective magnifier fits nicely in any handbag, shirt or trouser pocket and can also be used with a lanyard.

- Magnification from 3X – 14X (12 D – 56 D)
- Convenient integrated lanyard loop
- 1 magnifier with 6 different power variations

