

eMag 120/150 HD

Compact video magnifier in HD quality for mobile use!

True colours

With Color-Boost-Contrast



With unique
**Color-Boost-
Contrast**
for superior contrast
and brilliant colours!

 **SCHWEIZER**
Der LowVision-Spezialist



512002 eMag 120 HD

Magnification: 1.7x to 24x
 1.0x to 2.5x (self-view camera)
 Weight: 4.3 kg – with battery, 3.8 kg – without battery
 Dimensions (W x D x H): approx. 31.7 x 34.3 x 6.1 cm
 Screen size: 12 inch (30 cm)
 with 4-metre power cord and carrying case

515002 eMag 150 HD

Magnification: 2.3x to 31x
 1.25x to 3x (self-view camera)
 Weight: 4.8 kg – with battery, 4.3 kg – without battery
 Dimensions (W x D x H): approx. 40.6 x 34.3 x 6.1 cm
 Screen size: 15 inch (38 cm)
 with 4-metre power cord and carrying case



Front camera for self-viewing



Small and handy to store away



Ideal for desk work

Compact video magnifier

■ eMag 120/150 HD

Compact portable video magnifier in user-friendly design. Features all essential functions of a traditional desktop video magnifier. With unique Color-Boost-Contrast for superior contrast and brilliant colours – in HD quality! Coloured control elements ensure intuitive handling. Folds away easily for transportation.

- High-resolution HD auto-focus cameras for close-up view and self-view mode
- 2 viewing modes: document reading and self-view
- Individual magnification and colour settings are saved for each view
- Multiple comfort features including freeze frame, reading masks (black and shadow), reading lines
- 33 colour modes: true full colours, high-contrast full colours (Color-Boost-Contrast), black on white, white on black, yellow on blue, yellow on black and 27 additional colour modes
- Large controls in different colours facilitate intuitive handling
- Find feature
- Snapshot feature (images can be stored to SD card and played as slideshow)
- Space saving: folds quickly and easily for storage
- Monitor adjusts to any height or vertical tilt for optimal head posture and relaxed reading



■ Optional:

512012 Battery for eMag 120 HD and eMag 150 HD

Lithium-ion battery pack
 Charge time: approx. 3.5 h
 Operating time: eMag 120 HD approx. 4 h, eMag 150 HD approx. 3.5 h

510012 eMag Table

Trapezoid table with lockable castors, height-adjustable, dimensions (W x D x H): approx. 80 x 60 x 65-91 cm

512042 XY reading table for eMag 120 HD and eMag 150 HD

Dimensions: approx. 32 x 31.5 x 2.5 cm